

## NON-CONTACT FOREHEAD THERMOMETER

Forehead mode: 93.2 to 108°F/34 to 42.2°C

Surface mode: -7.6 to +176°F/-22 to +80°C

### TIPS

- If there is any temperature difference between the places where the device is stored and where it will be used, allow the thermometer to adjust to the room temperature where the subject is for at least 15 minutes before measurement.
- Remove hats or head bands before taking temperature.
- Do not take temperatures outside. If someone is just coming in from outside please wait at least 5 minutes to take temperature (outside temperature will affect the reading).
- **Unit must be within 1.5 inches or less** before taking temperature.
- Press and hold down the START button. When you hear one beep release the START button and you will have the temperature reading. When you hear two beeps (about 5 seconds) you can take the next temperature reading.
- NORMAL BODY TEMPERATURE IS BETWEEN 97 AND 99.5 °F.
- To silence alert, press and hold the POWER button (🔌) and release when the silent icon (🔇) flashes on the display.
- Gently clean the sensor lens periodically with wipe or cotton swab dampened with 70% or more isopropyl alcohol. Be gentle. Allow the lens to dry fully for at least 1 minute before using the thermometer.

