

Thermometer Probe Wipes – Packets

Model: PW200

Contains 200 individually packaged single-use wipes, saturated with 70% isopropyl alcohol to disinfect and sanitize hard, non-porous surfaces and help prevent cross contamination.



| | |
|-----------------------------|--|
| Features | 1.2 x 2.2 (in)/3 x 5.5 (cm) pads |
| Certification | FDA-listed, EPA & NQA certified |
| Antimicrobial Protection | Yes |
| Active Ingredient | Isopropanol (isopropyl alcohol) 70% |
| SARS-CoV-2 (COVID-19) | To kill SARS-CoV-2 (COVID-19), follow disinfection directions for the following virus(es): Rhinovirus |
| Disinfection & Sanitization | Tested in accordance with the AOAC Use Dilution test method [modified in the presence of 5% blood serum] |
| Bacteria Kill Time | 5 minutes: Staphylococcus aureus (ATCC 6538), Pseudomonas aeruginosa (ATCC 15442), Brevundimonas vesicularis (ATCC 51560), Burkholderia cepacia (ATCC 25416), Corynebacterium jeikeium (ATCC 43216), Escherichia coli (ATCC 11229), Micrococcus luteus (ATCC 49732), Staphylococcus epidermidis (ATCC 12228) 3 minutes: Staphylococcus aureus (ATCC 6538), Klebsiella pneumoniae (ATCC 4352), Enterobacter aerogenes (ATCC 13048) |
| Virucide | Tested in accordance with the Virucidal Qualification test method [modified in the presence of 5% blood serum] |
| Virus Kill Time | 5 minutes: Herpes Simplex Type I Virus (Herpes) (ATCC VR-733), Herpes Simplex Type II Virus (Genital Herpes) (ATCC VR-734), Influenza A Hong Kong H3N2 Virus (ATCC VR-544), Rhinovirus Type 37 (ATCC VR-1147) |
| Fungicide | Tested in accordance with the AOAC Fungicidal test method [modified in the presence of 5% blood serum] |
| Fungus Kill Time | 5 minutes: Candida albicans (ATCC 10231), Penicillium chrysogenum (ATCC 10003); 10 minutes: Trichophyton mentagrophytes (ATCC 9533) |
| Formulation Type | Ready-to-use |
| Surface Types | Hard non-porous (HN) |
| Use Sites | Healthcare, Institutional |
| Caution | Flammable; Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. |
| Weight | 7.3 oz / 207 g |
| Warranty | Expiration date on package in lieu of warranty |